

100 series Pump Heads only

Specification Sheet (last updated Jan 2010)



The 100 series, micro range of peristaltic pumps offer OEM's a low cost, compact and reliable means of dosing a wide range of liquids. They are available in various configurations of tubing type/ID bore as detailed below.

The easiest way to place an order or check prices is via our online shop at <http://www.williamson-shop.co.uk> but please contact us for any requirements not listed or to place an order over the phone.

Model No	Series	ml/rev*	Rollers	Tube id/material
100.PH.008/2	100 series pumphead	0.0300	2	0.8mm silicone
100.PH.008/4	100 series pumphead	0.0220	4	0.8mm silicone
100.PH.016/2	100 series pumphead	0.1000	2	1.6mm silicone
100.PH.016/4	100 series pumphead	0.0800	4	1.6mm silicone.
100.PH.030/2	100 series pumphead	0.2600	2	3.0mm silicone.
100.PH.030/4	100 series pumphead	0.2200	4	3.0mm silicone.
100D-008/2-KIT-4MM	100 double pumphead	0.0300	2	0.8mm silicone
100D-008/4-KIT-4MM	100 double pumphead	0.0220	4	0.8mm silicone
100D-016/2-KIT-4MM	100 double pumphead	0.1000	2	1.6mm silicone
100D-016/4-KIT-4MM	100 double pumphead	0.0800	4	1.6mm silicone.
100D-030/2-KIT-4MM	100 double pumphead	0.2600	2	3.0mm silicone.
100D-030/4-KIT-4MM	100 double pumphead	0.2200	4	3.0mm silicone.
101.PH.008/2	100 series pumphead	0.0275	2	0.8mm Norprene
101.PH.008/4	100 series pumphead	0.0200	4	0.8mm Norprene
101.PH.016/2	100 series pumphead	0.0900	2	1.6mm Norprene
101.PH.016/4	100 series pumphead	0.0700	4	1.6mm Norprene
101.PH.030/2	100 series pumphead	0.2000	2	3.0mm Norprene
101.PH.030/4	100 series pumphead	0.1600	4	3.0mm Norprene
101D-008/2-KIT-4MM	100 double pumphead	0.0275	2	0.8mm Norprene
101D-008/4-KIT-4MM	100 double pumphead	0.0200	4	0.8mm Norprene
101D-016/2-KIT-4MM	100 double pumphead	0.0900	2	1.6mm Norprene
101D-016/4-KIT-4MM	100 double pumphead	0.0700	4	1.6mm Norprene
101D-030/2-KIT-4MM	100 double pumphead	0.2000	2	3.0mm Norprene
101D-030/4-KIT-4MM	100 double pumphead	0.1600	4	3.0mm Norprene
102.PH.008/2	100 series pumphead	0.0250	2	0.8mm Viton
102.PH.008/4	100 series pumphead	0.0160	4	0.8mm Viton.
102.PH.016/2	100 series pumphead	0.0800	2	1.6mm Viton
102.PH.016/4	100 series pumphead	0.0600	4	1.6mm Viton.
102.PH.030/2	100 series pumphead	0.1800	2	3.0mm Viton
102.PH.030/4	100 series pumphead	0.1400	4	3.0mm Viton.
102D-008/2-KIT-4MM	100 double pumphead	0.0250	2	0.8mm Viton
102D-008/4-KIT-4MM	100 double pumphead	0.0160	4	0.8mm Viton.
102D-016/2-KIT-4MM	100 double pumphead	0.0800	2	1.6mm Viton
102D-016/4-KIT-4MM	100 double pumphead	0.0600	4	1.6mm Viton.
102D-030/2-KIT-4MM	100 double pumphead	0.1800	2	3.0mm Viton

*To calculate flow per minute, multiply ml/rev by RPM of motor gearbox being used. (on double pump heads ml/rev is per tube). Flow figures should be used only as an estimate. Customers should run their own tests to establish pump suitability

With models in blue, users should take care not to run the pump too fast. For advice please contact our technical department.

The Williamson Manufacturing Company Ltd provides no warranty on usage of pumps. We recommend that life tests be carried out prior to use. This information is given in good faith and believed to be correct and current at time of publishing. The Williamson Manufacturing Company Ltd cannot accept responsibility for its inaccuracy or any errors or omissions contained herein. Copyright The Williamson Manufacturing Company Ltd 2008